

User Guide

Relaxation Profile Protocol

This user guide has been created to educate and inform the reader about the Relaxation Profile Protocol.

For more information about NeXus, our BioTrace+ software, please visit our website or contact us.

www.mindmedia.com

The copyright of this document remains with Mind Media BV © 2020, and the contents of it may not be altered or copied. However we do permit unlimited distribution by electronic means of this document in the unaltered digital PDF format in which it is supplied. This document is not intended to replace scientific and clinical literature.

Contents

- Introduction..... 3**
- Required equipment 3**
- Measurement setup..... 4**
- Using the Relaxation Profile Protocol..... 5**
- Analysis and export to Relaxation Profile Template 10**
- Appendix: Installation..... 13**

Introduction

The Relaxation Profile User Guide provides a step-by-step review of how to install the Relaxation Profile, how to use the Relaxation Profile and how to exporting the data to the Relaxation Profile Template. This Relaxation Profile has been developed in cooperation with Inna Khazan, PhD, BCB.

Required equipment

Depending on the chosen setup, the following is required to perform the Stress test:

- Nexus-4, NeXus-10 or NeXus-32
- Skin Conductance Sensor
- Skin Conductance electrodes (Ag/AgCl)
- Temperature Sensor
- Blood Volume Pulse Sensor
- Respiration Sensor
- EXG Sensor
- EXG Ground
- Pre-gelled EMG electrodes*

*High quality electrodes like the Meditrace or ARBO electrodes are recommended to ensure good signal quality.

Measurement setup

Before the actual measurement can start, the equipment has to be connected. Detailed information on setting up the NeXus can be found in the NeXus User Manual or Quick Start.

Connect the sensors to the right NeXus inputs. Make sure the red dot of the connector is facing downward with the NeXus-10 or upward with the NeXus-32. Detailed information about sensor placement and preparation can be found in the measurement setup user guides.



Connect the EXG Ground to the Ground (Gnd) of the NeXus.

NeXus-10

- | | | |
|----------------|---|--------------------|
| C&D |  | EMG (EXG sensor) |
| E |  | Skin Conductance |
| F |  | Temperature |
| G |  | Blood Volume Pulse |
| H |  | Respiration |



NeXus-32




- | | | |
|------------------|---|--------------------|
| 27&28 |  | EMG (EXG sensor) |
| 29 |  | Skin Conductance |
| 30 |  | Blood Volume Pulse |
| 31 |  | Respiration |
| 32 |  | Temperature |

Using the Relaxation Profile Protocol

Make sure to install the Relaxation Profile Protocol first (Appendix: Installation).

The protocol will go through the following sequence.

Screen	Duration	Segment	Marker	
Baseline	120 seconds	Baseline	Baseline	
Relaxation Technique Breathing	180 seconds	Train	Breathing	
Pause	10 seconds		Pause	
Relaxation Technique Passive Muscle Relaxation	180 seconds	Train	Passive M. Relaxation	
Pause	10 seconds		Pause	


Relaxation Technique Progressive Muscle Relaxation	180 seconds	Train	Progressive M. Relaxation	
Pause	10 seconds		Pause	
Relaxation Technique Autogenic Training	180 seconds	Train	Autogenic Training	
Pause	10 seconds		Pause	
Relaxation Technique Guided Imagery	180 seconds	Train	Guided Imagery	
Finished			Finished	

Each segment will automatically be color coded in the session overview, including markers describing the relaxation technique for that particular segment.


Open BioTrace+.



Select 'Screen Browser'.




Select 'My Screens'.




Select 'Relaxation Protocol Inna Khazan'.

Double Click 'Relaxation Protocol' screen to open the therapist screen.



The therapist screen will be opened.



The following is displayed in the therapist screen: EMG (1); EMG (2); Temperature (3); Skin Conductance (4); Heart Rate (5); Breathing pattern (6); Blood Volume Pulse (7).

Turn on the NeXus device.

Start a recording.




Visually inspect the recorded signal to identify and prevent artifacts.

The actual Relaxation Profile protocol is started by clicking the 'Start protocol' button in the right top corner.




A dual monitor setup is recommended. Press the Windows logo key +P for extending display to dual monitor setup. The protocol will be opened automatically on the second monitor. When using a single monitor setup the protocol will be opened on the first screen.

Optionally use the 'Continue' button to proceed to the next screen in the protocol.



Once the protocol is finished, the recording will stop.




Click 'Yes' to save the session and to enter a description of the session. Confirm by clicking 'OK', the session is now saved.

Analysis and export to Relaxation Profile Template

Data can be reviewed in the *Session Overview* by clicking the following button in the Session Control Bar. Statistics of the *Session Overview* can be copied to the Relaxation Profile Template.



Right-click the *Session Overview* screen and choose 'Intra-session trend'.



The *Trend Report* screen will appear.

Select the 'Mean', 'Min. & Max.', 'Std.Dev.' and 'HRV' output options for all statistics.

Trend Report: Intra-Session (based on Segments)

This function computes an Intra-Session trend using segment-definitions. At least 2 segments must be defined. Choose your output options below. You can use Copy/Paste for trend data, with the editable text option.

	5:Sensor-E:SC/GSR	6:Sensor-F:Ti
Mean. value trend:		
1:Baseline	10.993	93.720
2:Train	11.245	94.105
3:Train	11.988	94.715
4:Train	11.426	94.109
5:Train	11.769	94.188
6:Train	10.906	94.339
Min. value trend:		
1:Baseline	10.258	93.113
2:Train	10.814	93.671
3:Train	10.619	94.210
4:Train	10.841	93.988
5:Train	10.683	94.048
6:Train	10.434	94.069

Mean.
 Min. & Max.
 Std.Dev.
 HRV (rms/sd)

OK 1 Text: editable text (copy/paste) Cancel

Select all data in the Trend Report by dragging the left mouse button.

Trend Report: Intra-Session (based on Segments)

This function computes an Intra-Session trend using segment-definitions. At least 2 segments must be defined. Choose your output options below. You can use Copy/Paste for trend data, with the editable text option.

	5:Sensor-E:SC/GSR	6:Sensor-F:Ti
Mean. value trend:		
1:Baseline	10.993	93.720
2:Train	11.245	94.105
3:Train	11.988	94.715
4:Train	11.426	94.109
5:Train	11.769	94.188
6:Train	10.906	94.339
Min. value trend:		
1:Baseline	10.258	93.113
2:Train	10.814	93.671
3:Train	10.619	94.210
4:Train	10.841	93.988
5:Train	10.683	94.048
6:Train	10.434	94.069

Mean.
 Min. & Max.
 Std.Dev.
 HRV (rms/sd)

OK 1 Text: editable text (copy/paste) Cancel

Trend Report: Intra-Session (based on Segments)

This function computes an Intra-Session trend using segment-definitions. At least 2 segments must be defined. Choose your output options below. You can use Copy/Paste for trend data, with the editable text option.


	5:Sensor-E:SC/GSR	6:Sensor-F:Ti
HRV RMSSD trend:		
1:Baseline	41.879	41.879
2:Train	28.826	28.826
3:Train	27.484	27.484
4:Train	47.210	47.210
5:Train	45.138	45.138
6:Train	40.360	40.360
HRV SDNN trend:		
1:Baseline	68.555	68.555
2:Train	44.496	44.496
3:Train	46.514	46.514
4:Train	78.161	78.161
5:Train	72.902	72.902
6:Train	66.502	66.502

Mean.
 Min. & Max.
 Std.Dev.
 HRV (rms/sd)

OK 1 Text: editable text (copy/paste) Cancel

Right click on the selected data and choose 'Copy'.


Open the Relaxation Profile Excel Template and right click on the first cell (A1) in the 'Data' tab.




Choose 'Paste' to copy all statistics to the Relaxation Profile Excel Template.

	A	B	C
1		5:Sensor-E:SC/GSR	6:Sensor-F:Temp.
2			
3	Mean. value trend:		
4	1:Baseline	10.993	93.720
5	2:Train	11.168	94.615
6	3:Train	12.246	94.346
7	4:Train	11.319	94.150
8	5:Train	11.274	94.242
9	6:Train	11.247	94.857
10			
11			
12	Min. value trend:		
13	1:Baseline	10.258	93.113
14	2:Train	10.619	93.905
15	3:Train	11.152	93.988
16	4:Train	10.683	93.993
17	5:Train	10.463	94.062
18	6:Train	10.428	94.628
19			
20			
21	Max. value trend:		
22	1:Baseline	10.660	94.240

Under the second tab 'Relaxation Profile Charts' tab, the mean statistics will be displayed in bar graphs.





Under the third tab 'Stats' the mean statistics will be displayed in tables.


Data | Relaxation Profile Charts | **Stats**

	Baseline 0-2min	Breathing 2-5min	Progressive Muscle Relaxation 5-8min	Passive Muscle Relaxation 8-11min	Autogenic Training 11-14min	Imagery 14-17min
HR						
Average	70,76	69,80	70,47	71,32	70,95	72,41
SDNN						
Average	93,11	44,50	52,32	82,49	72,90	67,72
Breathing Rate (BPM)						
Average	9,88	11,95	6,73	6,69	6,68	6,69
HF						
%HF	26,77	18,85	9,64	8,53	8,09	7,42
LF						
%LF	83,00	81,28	90,54	90,74	90,82	91,86
HRV Max-Min						
Average	9,58	11,22	17,89	16,72	15,10	17,04
Skin Conductance						
Average	10,99	11,17	12,25	11,32	11,27	11,25
Temperature						
Average	93,72	94,62	94,35	94,15	94,24	94,86
EMG A _ Right Trap						
Average	1,47	1,56	1,65	1,77	1,89	1,94
EMG B _ Left Trap						
Average	1,42	1,49	1,50	1,68	1,72	1,62


Appendix: Installation

Close the BioTrace+ software.

Open the '**NX Relaxation Profile.exe**'-file. The software will now search for BioTrace+ installations on the computer.



Select the BioTrace+ directory where the protocol needs to be installed (this step is only necessary when there are multiple versions of BioTrace+ installed) and click 'Next'.



Click 'Exit' once installation has finished.

