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EEG

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Appendix 4: Linked-ear cable setup

EXG sensor setup with disc electrodes

The following is required to perform EEG measurements in an EXG sensor setup with disc electrodes:

- Nexus-4, NeXus-10 or NeXus-32
- EXG Sensor
- EXG Ground
- EEG Discs* (Ag / AgCl)
- Nuprep
- Electrode paste (e.g. Ten20)
- Ear clips for EEG discs
- Measuring tape
- Tissue
- Cotton swabs
- Measuring tape
- Skin marker
- Scoop

*Recommended use is 10-20 times. The layer of Ag / AgCl will gradually wear off influencing signal quality.

One position EEG

Connect the EXG Sensor to input A&B of the NeXus-4 or NeXus-10 and input 25&26 of the NeXus-32. Make sure the red dot of the connector is facing downward with the NeXus-4 or NeXus-10 or upward with the NeXus-32.

Connect the EXG Ground to the Ground (Gnd) of the NeXus.

Sensors can be disconnected by pulling the silver ribbed part of the connector backward.

Place a disc electrode on to the red snap-on of pair one of the EXG sensor.

Place the black snap-on of pair one of the EXG sensor and the ground snap-on together with disc electrodes on to an ear clip.

Determine electrode positions according to the 10-20 system by using measuring tape and marking positions with a skin marker.

Rub the skin with a cotton swab and Nuprep at all electrode positions. Remove residual gel with a tissue. As an alternative to using Nuprep, alcohol pads are sometimes used.

Use the clip of the EXG cable to attach the EXG sensor to prevent pulling of the cables.

Apply some electrode paste (e.g. Ten20) to the electrodes by using for example a scoop.

Place the ear clip with the black (reference) snap-on with electrode on to the earlobe of the same side (ipsilateral) of the location of interest.

Place the earclip with ground snap-on with electrode on to the other earlobe.

Place the red (positive) snap-on with electrode on to the skin at the location of interest.

Minicap

Optionally use the Minicap or Microcap to fix the electrodes. The minicap also allows for placing the electrodes on the mastoids instead of the ear.

Now go to chapter 'BioTrace+' for further instructions.

Two positions EEG

Connect the EXG Sensor to input A&B of the NeXus-4 or NeXus-10 and input 25&26 of the NeXus-32. Make sure the red dot of the connector is facing downward with the NeXus-4 or NeXus-10 or upward with the NeXus-32.

To create one reference for both positive electrodes, use the linked ear cable.

In a linked-ear setup, the positive electrodes are referenced to two linked negative electrodes. In order to create a linked-ear setup, a Linked Ear cable is placed in between the EXG sensor and the input of the Nexus. Refer to the linked ear manual for an extensive overview.

Connect the EXG Ground to the Ground (Gnd) of the NeXus.

Sensors can be disconnected by pulling the silver ribbed part of the connector backward.

Place disc electrodes on to the two red snap-ons of the EXG sensor and place a disc electrode on to the EXG ground snap-on.

Place the two black snap-ons of the EXG sensor together with disc electrodes on to an ear clip.

Determine electrode positions according to the 10-20 system by using measuring tape and marking positions with a skin marker.

Rub the skin with a cotton swab and Nuprep at all electrode positions. Remove remaining gel with a tissue. As an alternative to using Nuprep, alcohol pads are sometimes used.

Use the clip of the EXG cable to attach the EXG sensor to prevent pulling of the cables.

Apply some electrode paste (e.g. Ten20) to the electrodes by using for example a scoop.

Place the ear clips with the black (reference) snap-on with electrode on to the left and right earlobe.

Place the ground electrode preferably near the other reference electrodes on the scalp.

Place the red (positive) snap-ons with electrode on to the skin at the location of interest.

Minicap

Optionally use the Minicap or Microcap to fix the electrodes.

Now go to chapter 'BioTrace+' for further instructions.

Four positions EEG

When measuring four positions EEG, a linked-ear setup is recommended. Connect the linked-ear cables to input A&B and C&D on the NeXus-10 and input 25&26 and 27&28 on the NeXus-32. Make sure the red dot of the connector is facing downward with the NeXus-10 or upward with the NeXus-32.


Sensors can be disconnected by pulling the silver ribbed part of the connector backward.


Connect the EXG Ground to the Ground (Gnd) of the NeXus.

Connect the black connectors of the two Linked ear cables by using the small black tube.

Connect the EXG Sensors to the Linked Ear cables. Make sure the red dots of the connectors are aligned.

Place disc electrodes on to the four red snap-ons of the EXG sensor and place a disc electrode on to the EXG ground snap-on.

Place the two black snap-ons of the EXG sensor which is connected to input A&B of the NeXus-10 and input 25&26 of the NeXus-32 together with disc electrodes on to an ear clip.

Determine electrode positions according to the 10-20 system by using measuring tape and marking positions with a skin marker.

Rub the skin with a cotton swab and Nuprep at all electrode positions. Remove remaining gel with a tissue. As an alternative to using Nuprep, alcohol pads are sometimes used.

Use the clip of the EXG cable to attach the EXG sensor to prevent pulling of the cables.

Apply some electrode paste (e.g. Ten20) to the electrodes by using for example a scoop.

Place the ear clips with the black (reference) snap-on with electrode on to the left and right earlobe.

Place the ground electrode preferably near the other reference electrodes on the scalp.

Place the red (positive) snap-ons with electrode on to the skin at the location of interest.

Minicap

Optionally use the Minicap or Microcap to fix the electrodes.

Now go to chapter 'BioTrace+' for further instructions.

EXG sensor TP setup with Minicap electrodes

The following is required to perform EEG measurements in an EXG sensor TP setup with minicap electrodes:

- Nexus-4, NeXus-10 or NeXus-32
- EXG Sensor TP
- EXG Ground TP
- EEG Minicap Electrodes (sintered Ag/AgCl)
- Electro gel
- Ear clips for EEG Minicap electrodes
- Syringe with blunt needle
- Minicap or Microcap
- Tissue
- Cotton swabs
- Measuring tape
- Skin marker

One position EEG

Connect the EXG Sensor TP to input A&B of the NeXus-4 or NeXus-10 and input 25&26 of the NeXus-32. Make sure the red dot of the connector is facing downward with the NeXus-4 or NeXus-10 or upward with the NeXus-32.

Connect the EXG Ground TP to the Ground (Gnd) of the NeXus.


Sensors can be disconnected by pulling the silver ribbed part of the connector backward.

Place Minicap electrodes on to pair one of the EXG sensor TP and place a Minicap electrode on to the EXG ground TP. Slightly twist the connector in to EXG TP connection in order to make good connection.

Place the black TP connection (reference) with electrode of pair one and the ground with electrode on to an ear clip.

Determine electrode positions according to the 10-20 system by using measuring tape and marking positions with a skin marker.

Use the clip of the EXG TP cable to attach the EXG TP sensor to prevent pulling of the cables.

Place the Minicap or Microcap on the head.

Apply some Electro gel to the electrodes by using a cotton swab or a syringe with blunt needle.

Due to the fluid nature of electro gel, rubbing the skin with an abrasive gel is usually not necessary.

When using electrode paste (e.g. Ten20), rub the skin with a cotton swab and Nuprep at all electrode positions. The location of the ground electrode is preferably near other electrodes on the scalp. Remove residual gel with a tissue. As an alternative to using Nuprep, alcohol pads are sometimes used.

Place the ear clip with the black (reference) electrode on to the earlobe of the same side (ipsilateral) of the location of interest.

Place the earclip with ground snap-on with electrode on to the other earlobe.

Place the red (positive) electrode on to the skin at the location of interest. Use the Minicap or Microcap to fix the electrode. The minicap also allows for placing the electrodes on the mastoids instead of the ear. In that case the earclips aren't necessary.

Now go to chapter 'BioTrace+' for further instructions.

Two positions EEG

Connect the EXG Sensor TP to input A&B of the NeXus-4 or NeXus-10 and input 25&26 of the NeXus-32. Make sure the red dot of the connector is facing downward (NeXus-4 or NeXus-10) or upward (NeXus-32).

To create one reference for both positive electrodes, use the linked ear cable.

In a linked-ear setup, the positive electrodes are referenced to two linked negative electrodes. In order to create a linked-ear setup, a Linked Ear cable is placed in between the EXG sensor and the input of the Nexus. Refer to the linked ear manual for an extensive overview.

Connect the EXG Ground TP to the Ground (Gnd) of the NeXus.


Sensors can be disconnected by pulling the silver ribbed part of the connector backward.

Place Minicap electrodes on the EXG sensors TP and place a Minicap electrode on to the EXG ground TP. Slightly twist the connector in to EXG TP connection in order to make good connection.

Place the black TP connection (reference) electrodes on to an ear clip.

Determine electrode positions according to the 10-20 system by using measuring tape and marking positions with a skin marker.

Use the clip of the EXG TP cable to attach the EXG TP sensor to prevent pulling of the cables.

Place the Minicap or Microcap on the head.

Apply some Electro gel to the electrodes by using a cotton swab or a syringe with blunt needle.

Due to the fluid nature of electro gel, rubbing the skin with an abrasive gel is usually not necessary.

When using electrode paste (e.g. Ten20), rub the skin with a cotton swab and Nuprep at all electrode positions. The location of the ground electrode is preferably near other electrodes on the scalp. Remove residual gel with a tissue. As an alternative to using Nuprep, alcohol pads are sometimes used.

Place the ear clips with the black (reference) electrodes on to the left and right earlobe.

Place the ground electrode preferably near the other reference electrodes on the scalp. Use the Minicap or Microcap to fix the electrodes.

Place the red (positive) electrodes on to the skin at the location of interest. Use the Minicap or Microcap to fix the electrode. The minicap also allows for placing the electrodes on the mastoids instead of the ear. In that case the earclips aren't necessary.

Now go to chapter 'BioTrace+' for further instructions.

Four positions EEG

When measuring four positions EEG, a linked-ear setup is recommended. Connect the linked-ear cables to input A&B and C&D on the NeXus-10 and input 25&26 and 27&28 on the NeXus-32. Make sure the red dot of the connector is facing downward (NeXus-10) or upward (NeXus-32).


Sensors can be disconnected by pulling the silver ribbed part of the connector backward.

Connect the EXG Ground TP to the Ground (Gnd) of the NeXus.

Connect the black connectors of the two Linked ear cables by using the small black tube.


Connect the EXG Sensors TP to the Linked Ear cables. Make sure the red dots of the connectors are aligned.

Place Minicap electrodes on to the red EXG sensors TP and place a Minicap electrode on to the EXG ground TP. Slightly twist the connector in to EXG TP connection in order to make good connection.

Place black Minicap electrodes on to the black electrodes of the EXG sensor TP which is connected to input A&B of the NeXus-10 and input 25&26 of the NeXus-32.

Place the black TP connection (reference) electrodes on to an ear clip.

Determine electrode positions according to the 10-20 system by using measuring tape and marking positions with a skin marker.

Use the clip of the EXG TP cable to attach the EXG TP sensor to prevent pulling of the cables.

Place the Minicap or Microcap on the head.

Apply some Electro gel to the electrodes by using a cotton swab or a syringe with blunt needle.

Due to the fluid nature of electro gel, rubbing the skin with an abrasive gel is usually not necessary.

When using electrode paste (e.g. Ten20), rub the skin with a cotton swab and Nuprep at all electrode positions. The location of the ground electrode is preferably near other electrodes on the scalp. Remove residual gel with a tissue. As an alternative to using Nuprep, alcohol pads are sometimes used.

Place the ear clips with the black (reference) electrodes on to the left and right earlobe.

Place the ground electrode preferably near the other reference electrodes on the scalp. Use the Minicap or Microcap to fix the electrodes.

Place the red (positive) electrodes on to the skin at the location of interest. Use the Minicap or Microcap to fix the electrode. The minicap also allows for placing the electrodes on the mastoids instead of the ear. In that case the earclips aren't necessary.

Now go to chapter 'BioTrace+' for further instructions.

EXG sensor TP setup with cup or ring electrodes

The following is required to perform EEG measurements in an EXG sensor TP setup with cup or ring electrodes:

- Nexus-4, NeXus-10 or NeXus-32
- EXG Sensor TP
- EXG Ground TP
- EEG cup or ring electrodes (Ag/AgCl)
- Nuprep
- Electrode paste (e.g. Ten20)
- Measuring tape
- Tissue
- Cotton swabs
- Skin marker
- Scoop

One position EEG

Connect the EXG Sensor TP to input A&B of the NeXus-4 or NeXus-10 and input 25&26 of the NeXus-32. Make sure the red dot of the connector is facing downward with the NeXus-4 or NeXus-10 or upward with the NeXus-32.


Connect the EXG Ground TP to the Ground (Gnd) of the NeXus.

Sensors can be disconnected by pulling the silver ribbed part of the connector backward.

Place cup or ring electrodes on to pair one of the EXG sensor TP and place a cup or ring electrode on to the EXG ground TP.

Determine electrode positions according to the 10-20 system by using measuring tape and marking positions with a skin marker.

Rub the skin with a cotton swab and Nuprep at all electrode positions. Remove residual gel with a tissue. As an alternative to using Nuprep, alcohol pads are sometimes used.

Use the clip of the EXG TP cable to attach the EXG TP sensor to prevent pulling of the cables.

Apply some electrode paste (e.g. Ten20) to the electrodes by using for example a scoop.

Place the black (reference) electrode on to the earlobe of the same side (ipsilateral) of the location of interest.

Place the ground electrode on to the other earlobe.

Place the red (positive) electrode on to the skin at the location of interest.

Now go to chapter 'BioTrace+' for further instructions.

Two positions EEG

Connect the EXG Sensor TP to input A&B of the NeXus-4 or NeXus-10 and input 25&26 of the NeXus-32. Make sure the red dot of the connector is facing downward (NeXus-4 or NeXus-10) or upward (NeXus-32).

To create one reference for both positive electrodes, use the linked ear cable.

In a linked-ear setup, the positive electrodes are referenced to two linked negative electrodes. In order to create a linked-ear setup, a Linked Ear cable is placed in between the EXG sensor and the input of the Nexus. Refer to the linked ear manual for an extensive overview.

Connect the EXG Ground TP to the Ground (Gnd) of the NeXus.

Sensors can be disconnected by pulling the silver ribbed part of the connector backward.

Place cup or ring electrodes on the EXG sensors TP and place a cup or ring electrode on to the EXG ground TP.

Determine electrode positions according to the 10-20 system by using measuring tape and marking positions with a skin marker.

Rub the skin with a cotton swab and Nuprep at all electrode positions. Remove remaining gel with a tissue. As an alternative to using Nuprep, alcohol pads are sometimes used.

Use the clip of the EXG TP cable to attach the EXG TP sensor to prevent pulling of the cables.

Apply some electrode paste (e.g. Ten20) to the electrodes by using for example a scoop.

Place the black (reference) electrodes on to the left and right earlobe.

Place the ground electrode preferably near the other reference electrodes on the scalp.

Place the red (positive) electrodes on to the skin at the locations of interest.

Now go to chapter 'BioTrace+' for further instructions.

Four positions EEG

When measuring four positions EEG, a linked-ear setup is recommended. Connect the linked-ear cables to input A&B and C&D on the NeXus-10 and input 25&26 and 27&28 on the NeXus-32. Make sure the red dot of the connector is facing downward (NeXus-10) or upward (NeXus-32).


Sensors can be disconnected by pulling the silver ribbed part of the connector backward.

Connect the EXG Ground TP to the Ground (Gnd) of the NeXus.

Connect the black connectors of the two Linked ear cables by using the small black tube.

Connect the EXG Sensors TP to the Linked Ear cables. Make sure the red dots of the connectors are aligned.

Place cup or ring electrodes on to the red EXG sensors TP and place a cup or ring electrode on to the EXG ground TP.

Place cup or ring electrodes on to the black electrodes of the EXG sensor TP which is connected to input A&B of the NeXus-10 and input 25&26 of the NeXus-32.

Determine electrode positions according to the 10-20 system by using measuring tape and marking positions with a skin marker.

Rub the skin with a cotton swab and Nuprep at all electrode positions. Remove remaining gel with a tissue. As an alternative to using Nuprep, alcohol pads are sometimes used.

Use the clip of the EXG TP cable to attach the EXG TP sensor to prevent pulling of the cables.

Apply some electrode paste (e.g. Ten20) to the electrodes by using for example a scoop.

Place the black (reference) electrodes on to the left and right earlobe.

Place the ground electrode preferably near the other reference electrodes on the scalp.

Place the red (positive) electrodes on to the skin at the locations of interest.

Now go to chapter 'BioTrace+' for further instructions.

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Das Urheberrecht an diesem Dokument, bleibt bei Mind Media BV © 2019 und der Inhalt darf nicht verändert oder kopiert werden. Wir erlauben jedoch die unbegrenzte elektronische Verbreitung dieses Dokuments im unveränderten digitalen PDF-Format, in dem es geliefert wird. Dieses Dokument ersetzt nicht die wissenschaftliche und klinische Literatur.